

Steam Sterilization (Autoclave Statement)

1. AMI Reusable Devices are compatible with all steam sterilization processes including Flash, Gravity and high Vac. Follow the sterilizer manufacturer's instructions for the various processes.
2. AMI LEEP/LLETZ coated instruments can be sterilized as described in the Operational Guidelines and Instructions For Use available in the IFU section.

Please note that although these instrument's high temperature Nylon coating is compatible with Flash Sterilization up to 275° F, we recommend to use 250°F sterilization process to increase the life of the coating on the instruments. ETO-Gas Sterilization, if available, will provide maximum life for the coating and is the recommended process.